



LABINCO L202 Rock Roller Variable Speed

Description

Additional Information

Number of positions	13
Max Quantity	1 kg
Motor Rating input(W)	30W
Motor Rating output(W)	25 W
Speed display	analog
Speed range(rpm)	1-25 rpm
Speed control	analog
Dimensions (W x D x H)	320 x 310 x 190 mm
Weight(kg)	3.6 kg (+ 1.2 kg control unit)
Permissible ambient temperature	5 - 50 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	50 W