




## LABINCO L209350 Rolling Bench 350 mm long rolls with strong motor

### Description

#### Roller bench L209350

- Analog unit with timer.
- Supplied as standard with 2 rolls.
- Expandable with 3 rolls to a max of 5
- 1 fixed role and 1 that you can place in multiple positions

*multi-functional; for multiple applications, from test tubes to large bottles*

- Rolling Bench to keep samples in motion.
- One roll is driven. Standard with 2 rolls, 1 fixed and 1 loose.
- roll length 350 mm. 46 mm. To order more rolls, Art.nr.35209
- You can easily place maximal 5 rolls on the rolling bench
- With timer and strong motor
- Speed 1-240 rpm
- Heart to heart distance from 7 mm up to 2150 mm
- Dimensions basic unit (WxDxH) 580x 350 x 130 mm
- Weight 8.4 kg
- Art.nr.20935 (230V)

### Additional Information

|                       |           |
|-----------------------|-----------|
| Number of positions   | 1 up to 4 |
| Max Quantity          | 15 kg     |
| Motor Rating input(W) | 100W      |

|                                                   |                               |
|---------------------------------------------------|-------------------------------|
| <b>Motor Rating output(W)</b>                     | 75 W                          |
| <b>Speed display</b>                              | analog                        |
| <b>Speed range(rpm)</b>                           | 1-240 rpm                     |
| <b>Stirring bar length max(mm)</b>                | digital 1s. up to 99h 59m 59s |
| <b>Speed control</b>                              | analog                        |
| <b>Dimensions (W x D x H)</b>                     | 580 x 350 x 130 mm            |
| <b>Weight(kg)</b>                                 | 8.4 kg                        |
| <b>Permissible ambient temperature</b>            | 5 - 50 C                      |
| <b>Permissible relative humidity</b>              | 80 %                          |
| <b>Protection class according to DIN EN 60529</b> | IP 42                         |
| <b>Voltage</b>                                    | 220 - 230 V                   |
| <b>Frequency</b>                                  | 50/60 Hz                      |
| <b>Power input</b>                                | 150W                          |