



LABINCO L209450 Rolling Bench for tanning salmon skins

Description

Additional Information

Number of positions	1 up to 4
Motor Rating input(W)	100 W
Motor Rating output(W)	75 W
Speed display	Analog
Speed range(rpm)	1 - 240 RPM
Speed control	Analog
Dimensions (W x D x H)	680 x350 x 130 mm
Weight(kg)	9.4 kg
Permissible ambient temperature	5 - 50 C
Permissible relative humidity	80%
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	150 W