




LABINCO L209650 Rolling Bench 650 mm long rolls with strong motor

Description

Roller bench L209650

- Analog unit with timer.
- Supplied as standard with 2 rolls.
- Expandable with 3 rolls to a max of 5
- 1 fixed role and 1 that you can place in multiple positions

multi-functional; for multiple applications, from test tubes to large bottles

- Rolling Bench to keep samples in motion.
- One roll is driven. Standard with 2 rolls, 1 fixed and 1 loose.
- roll length 650 mm. 46 mm. To order more rolls, Art.nr.65209
- You can easily place maximal 5 rolls on the rolling bench
- With timer and strong motor
- Speed 1-240 rpm
- Heart to heart distance from 7 mm up to 2150 mm
- Dimensions basic unit (WxDxH) 880x 350 x 130 mm
- Weight 11.5 kg
- Art.nr.20965 (230V)

Additional Information

Number of positions	1 up to 4
Max Quantity	15 kg
Motor Rating input(W)	100W

Motor Rating output(W)	75 W
Speed display	analog
Speed range(rpm)	1-240 rpm
Stirring bar length max(mm)	digital 1s. up to 99h 59m 59s
Speed control	analog
Dimensions (W x D x H)	880 x 350 x 130 mm
Weight(kg)	11.5 kg
Permissible ambient temperature	5 - 50 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	150W