



## LABINCO L22 Lum Magnetic Stirrer 5 ltr.

### Description

- Topplate 160 mm diam. Opal polycarbonate 10-1000 rpm
  - Stirring capacity 5 Ltr
  - LED-lighting
  - Dimensions (WxDxH) 200 x 250 x 85 mm
  - Weight 2.9 kg
  - Art.nr.22000 (230V) 22010 (115V)
  - Incl. spinbar
- 
- Topplate 160 mm diam. Opal polycarbonate 10-1000 rpm
  - Stirring capacity 5 Ltr
  - LED-lighting
  - Dimensions (WxDxH) 200 x 250 x 85 mm
  - Weight 2.9 kg
  - Art.nr.22000 (230V) 22010 (115V)
  - Incl. spinbar

### Additional Information

|                               |              |
|-------------------------------|--------------|
| <b>Number of positions</b>    | 1            |
| <b>Max Quantity</b>           | 5 ltr        |
| <b>Motor Rating input(W)</b>  | 50W          |
| <b>Motor Rating output(W)</b> | 20 W         |
| <b>Speed display</b>          | analog       |
| <b>Speed range(rpm)</b>       | 10 - 1000rpm |

|   |                        |
|---|------------------------|
| <b>Stirring bar length max(mm)</b>                | 60 mm                  |
| <b>Speed control</b>                              | stepless               |
| <b>Dimensions (W x D x H)</b>                     | 200 x 250 x 85 mm      |
| <b>Weight(kg)</b>                                 | 2.9 kg                 |
| <b>Permissible ambient temperature</b>            | 5 - 40 C               |
| <b>Permissible relative humidity</b>              | 80 %                   |
| <b>Protection class according to DIN EN 60529</b> | IP 42                  |
| <b>Voltage</b>                                    | 220 - 230 V (110/115V) |
| <b>Frequency</b>                                  | 50/60 Hz               |
| <b>Power input</b>                                | 75 W                   |