




LABINCO L26 Small Rotary Mixer

Description

Rotary-mixer Model L26

- Provides serial rotation of test-tube with diameter-range: 8 to 27 mm
- All aluminum construction
- Speed: stepless controllable: 5 to 50 rpm
- Quick coupling / exchange of rotating drum / bar
- Provided with stainless-steel drip-tray
- Small Rotary Mixer for serial rotation of test-tubes or drums
- Speed 1-50 rpm, analog
- Provided with stainless steel drip tray
- Quick coupling for exchange of rotating drum or bar
- Dimensions (WxDxH) 540 x 160 x 185 mm
- Weight 7.5 kg
- Art.nr.26000 (230V) 26010 (115V)
- See also the other models for more options on www.labinco-bv.com

Additional Information

Number of positions	1 bar or drum holder
Max Quantity	8 up to 27 mm or max 2 kg
Motor Rating input(W)	40 W
Motor Rating output(W)	25 W
Speed display	analog
Speed range(rpm)	1 - 50 rpm

Speed control	stepless
Dimensions (W x D x H)	540 x 160 x 185 mm
Weight(kg)	7.5 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	60 W