




LABINCO L28 Test-Tube Rotator

Description

Model ?? L28

- For mixing, diffusion, dialysing, extraction etc.
 - SPEED 10 rpm With strong, geared motor
 - Angle of rotation can be set between 90 and 27
 - Tubes may be exchanged while rotator is running
 - Comes with line-cord, but without disks
-
- Test-Tube Rotator for serial rotation of test-tubes 4-40 mm
 - Speed 10 rpm not adjustable
 - Angle of rotation can be set between 90 and 27
 - Dimensions (WxDxH) 200x 200 x 220 mm
 - Weight 2.6 kg
 - Art.nr.28000 (230V) 28010 (115V)
 - See also the other models for more options on www.labinco-bv.com

Additional Information

Number of positions	1 for max 2 disks
Max Quantity	1kg
Motor Rating input(W)	30W
Motor Rating output(W)	25W
Speed range(rpm)	10rpm
Speed control	not adjustable

Dimensions (W x D x H)	200 x 200 x 220 mm
Weight(kg)	2.6 kg
Permissible ambient temperature	5 - 60 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	35 W