



LABINCO L46 Power Mixer

Description

Additional Information

Number of positions	1
Max Quantity	Test tubes 7 mm up to 35 mm Dia.
Motor Rating input(W)	40 W
Motor Rating output(W)	25 W
Speed display	Analog
Speed range(rpm)	50 - 2500 rpm
Speed control	Stepless
Dimensions (W x D x H)	135 x 120 135 mm
Weight(kg)	3.2 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80%
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	50 W