



LABINCO L714 4 position Magnetic Stirrer 4 x 5 ltr.

Description

- Topplate 545 x 205 mm diam. Stainles steel 50-1250 rpm
 - 4 positions in one row
 - Stirring capacity 4 x 5 Ltr
 - Dimensions (WxDxH) 550 x 245 x 105 mm
 - Weight 8 kg
 - Art.nr.71400 (230V)
 - Incl. 4 x spinbar
-
- Topplate 545 x 205 mm diam. Stainles steel 50-1250 rpm
 - 4 positions in one row
 - Stirring capacity 4 x 5 Ltr
 - Dimensions (WxDxH) 550 x 245 x 105 mm
 - Weight 8 kg
 - Art.nr.71400 (230V)
 - Incl. 4 x spinbar

Additional Information

Number of positions	4 positions in one row
Max Quantity	4 x 5 ltr
Motor Rating input(W)	4 x 35W
Motor Rating output(W)	4 x 15 W
Speed display	analog
Speed range(rpm)	50-2500rpm

Stirring bar length max(mm)	60 mm
Speed control	stepless
Dimensions (W x D x H)	550 x 245 x 105 mm
Weight(kg)	8 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	150 W