




LABINCO LD208350 Rolling Bench 350 mm long rolls

Description

- With very strong motor. Up to 60 Kg
 - Touchscreen
 - Speedcontrol, up to 600 rpm
 - Timer settings
 - Voltage requirements 230 Volts 50/60 Hz.
-
- With very strong motor. Up to 60 Kg
 - Touchscreen
 - Speedcontrol, up to 600 rpm
 - Timer settings
 - Voltage requirements 230 Volts 50/60 Hz.

Additional Information

Motor Rating input(W)	250W
Motor Rating output(W)	175W
Speed display	Digital
Speed range(rpm)	1-600 rpm
Speed control	Digital
Dimensions (W x D x H)	580 x 350 x 130 mm
Weight(kg)	8,4 kg
Permissible ambient temperature	5 - 50 C
Permissible relative humidity	80%
Protection class according to DIN EN 60529	IP 42

No image found

Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	300W