




## LABINCO LD31 Large Hotplate 350°C

### Description

- Hotplate 460 x 350 mm Aluminium
- Hotplate temperature 30-350C digital
- Dimensions (WxDxH) 470 x 370 x 145 mm
- Weight 9.1 kg
- Temperature sensor on request
- Art.nr.31001 (230V)

- Hotplate 460 x 350 mm Aluminium
- Hotplate temperature 30-350C digital
- Dimensions (WxDxH) 470 x 370 x 145 mm
- Weight 9.1 kg
- Temperature sensor on request
- Art.nr.31001 (230V)

### Additional Information

|                                               |                                |
|-----------------------------------------------|--------------------------------|
| <b>Number of positions</b>                    | 1 large hotplate for more pos. |
| <b>Max Quantity</b>                           | 1 ml -50 ltr or more           |
| <b>Heat output</b>                            | 2000 W                         |
| <b>Heating rate</b>                           | 7 K/min                        |
| <b>Heating temperature range</b>              | room temp. - 350 C             |
| <b>Heat control</b>                           | digital                        |
| <b>Heat control accuracy</b>                  | 1 K                            |
| <b>Connection for ext. temperature sensor</b> | on request                     |

|                                            |                    |
|--------------------------------------------|--------------------|
| Control accuracy with sensor               | 1 K                |
| Temperatureconstancyinmedium               | 1 K                |
| Dimensions (W x D x H)                     | 470 x 370 x 145 mm |
| Weight(kg)                                 | 9.1 kg             |
| Permissible ambient temperature            | 5 - 40 C           |
| Permissible relative humidity              | 80 %               |
| Protection class according to DIN EN 60529 | IP 42              |
| Voltage                                    | 220 - 230 V        |
| Frequency                                  | 50/60 Hz           |
| Power input                                | 2000 W             |