



LABINCO Special big rolling bench

Description

Additional Information

| | |
|---|-----------------|
| Number of positions | 1 up to 8 |
| Motor Rating input(W) | 100 W |
| Motor Rating output(W) | 75 W |
| Speed display | Analog |
| Speed range(rpm) | 5 - 240 rpm |
| Speed control | Stepless |
| Dimensions (W x D x H) | 950 x 550 x 135 |
| Weight(kg) | 15.5 |
| Permissible ambient temperature | 5 - 40 C |
| Permissible relative humidity | 80% |
| Protection class according to DIN EN 60529 | IP 42 |
| Voltage | 220 - 230 V |
| Frequency | 50/60 Hz |
| Power input | 150 W |