



## LABINCO Stacked multi position magnet stirrer

### Description

This heating system are special made for a customer in Afrika. It is used to heat monsters in ceramic tubes with 6 different dimensions rapidly, up to 500 degrees.

All our custom made products are designed / developed and created in our own factory in The Netherlands.

Custom made product, let us know your wishes.

### Additional Information

<b>Number of positions</b>	12
<b>Speed range(rpm)</b>	1 - 25
<b>Voltage</b>	100 - 230 V
<b>Power input</b>	100 W