



## LABINCO very strong Power Mixer

### Description

### Additional Information

<b>Number of positions</b>	1
<b>Motor Rating input(W)</b>	90 W
<b>Motor Rating output(W)</b>	60 W
<b>Speed display</b>	Analog
<b>Speed range(rpm)</b>	50 - 2500 rpm
<b>Stirring bar length max(mm)</b>	Ø 45 mm or more (on request)
<b>Speed control</b>	Analog (stepless)
<b>Dimensions (W x D x H)</b>	200 x 270 x 170 mm
<b>Weight(kg)</b>	14 kg
<b>Permissible ambient temperature</b>	5 - 40 °C
<b>Permissible relative humidity</b>	80%
<b>Protection class according to DIN EN 60529</b>	IP 42
<b>Voltage</b>	100 - 230 V
<b>Frequency</b>	50/60 Hz
<b>Power input</b>	125 W