



LABINCO Flocculator 8 positions

Beschreibung

Custom made product, let us know your wishes.

Zustzliche Informationen

Number of positions	8
Max Quantity	10 ltr
Motor Rating input(W)	35 W
Motor Rating output(W)	25 W
Speed display	Analog
Speed range(rpm)	5 - 300 rpm
Stirring bar length max(mm)	80 mm
Speed control	Stepless (each position)
Dimensions (W x D x H)	1420 x 330 x 730 mm
Weight(kg)	48 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80%
Protection class according to DIN EN 60529	IP 42
Voltage	100 - 230 V
Frequency	50/60 Hz
Power input	250 W