




LABINCO Fresh cow's milk stirrer in refrigerator coffee machines

Beschreibung

This is a magnetic stirrer built into a refrigerator. This is to keep fresh milk moving at a low speed so the cream doesn't float. It is a separate refrigerator that sits next to a coffee machine and so you always have nice fresh milk with your cappuccino or latte.

All our custom made products are designed / developed and created in our own factory in The Netherlands.

Custom made product, let us know your wishes.

Zustätzliche Informationen

Number of positions	1 (or more on request)
Max Quantity	6 ltr.
Motor Rating input(W)	35 W
Motor Rating output(W)	25 W
Speed display	Analog
Speed range(rpm)	45 rpm
Stirring bar length max(mm)	80 mm
Speed control	On / Off
Dimensions (W x D x H)	On request
Weight(kg)	1 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80%
Protection class according to DIN EN 60529	IP 42

No image found

Voltage	100 - 230 V
Frequency	50/60 Hz
Power input	50 W