




LABINCO L17 Large Magnetic Stirrer 60 ltr.

Beschreibung

- Optical speed-detection
- Reversible stirring
- Topplate st.st.
- Micro-processor controlled
- Topplate 320 x 340 mm Stainless steel 50-1000 rpm
- Stirring capacity 60 Ltr
- Dimensions (WxDxH) 320 x 340 x 105 mm
- Art.nr.17000 (230V) 17010 (115V)
- Weight 5.2 kg
- Incl. spinbar

Zustzliche Informationen

Number of positions	1
Max Quantity	60 ltr
Motor Rating input(W)	50W
Motor Rating output(W)	20 W
Speed display	analog
Speed range(rpm)	50 - 1000rpm
Stirring bar length max(mm)	120 mm
Speed control	stepless
Dimensions (W x D x H)	320 x 340 x 105 mm
Weight(kg)	5.2 kg

Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115/110V)
Frequency	50/60 Hz
Power input	50 W