



LABINCO L22 Lum Magnetic Stirrer 5 ltr.

Beschreibung

- Topplate 160 mm diam. Opal polycarbonate 10-1000 rpm
 - Stirring capacity 5 Ltr
 - LED-lighting
 - Dimensions (WxDxH) 200 x 250 x 85 mm
 - Weight 2.9 kg
 - Art.nr.22000 (230V) 22010 (115V)
 - Incl. spinbar
-
- Topplate 160 mm diam. Opal polycarbonate 10-1000 rpm
 - Stirring capacity 5 Ltr
 - LED-lighting
 - Dimensions (WxDxH) 200 x 250 x 85 mm
 - Weight 2.9 kg
 - Art.nr.22000 (230V) 22010 (115V)
 - Incl. spinbar

Zustzliche Informationen

Number of positions	1
Max Quantity	5 ltr
Motor Rating input(W)	50W
Motor Rating output(W)	20 W
Speed display	analog
Speed range(rpm)	10 - 1000rpm

Stirring bar length max(mm)	60 mm
Speed control	stepless
Dimensions (W x D x H)	200 x 250 x 85 mm
Weight(kg)	2.9 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (110/115V)
Frequency	50/60 Hz
Power input	75 W