



LABINCO L224 4 position Lum Magnetic Stirrer 4 x 5 ltr.

Beschreibung

Model L224 LUM-Stirrer

- Glass topplate
- Electronic speed control
- With LED lighting
- Topplate 185 x 530 mm diam. Opal glass 10-1000 rpm
- Stirring capacity 4 x 5 Ltr
- LED-lighting
- Dimensions (WxDxH) 550 x 235 x 105 mm
- Weight 9 kg
- Art.nr.22400 (230V)
- Incl. 4 x spinbar

Zustzliche Informationen

Number of positions	4
Max Quantity	4 x 5 ltr
Motor Rating input(W)	4 x 35W
Motor Rating output(W)	4 x 15 W
Speed display	analog
Speed range(rpm)	10 - 1000rpm
Stirring bar length max(mm)	60 mm
Speed control	stepless

Dimensions (W x D x H)	550 x 235 x 105 mm
Weight(kg)	9 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	175 W