



LABINCO L23 Magnetic Stirrer 20 ltr.

Beschreibung

Model L23 Magnetic-Stirrer

- Electronic speed control
- Speed; 50 - 1250 rpm
- Housing: 2 piece casted aluminum
- Dripwater-spillage proved design and construction
- Electrostatically applied powder coated
- M-10 threaded support rod connection
- Supplied with te on magnetic stirring bar line cord and manual
- 10 liter capacity stainless steel topplate: 180 mm
- Topplate 175 mm diam. Stainless steel
- Stirring capacity 20 Ltr
- Dimensions (WxDxH) 200 x 250 x 90 mm
- Weight 2.9 kg
- Art.nr.23000 (230V) 23010 (115V)
- Incl. spinbar

Zustzliche Informationen

Number of positions	1
Max Quantity	20 ltr
Motor Rating input(W)	35 W
Motor Rating output(W)	15 W
Speed display	analog

Speed range(rpm)	50-2500rpm
Stirring bar length max(mm)	80 mm
Speed control	stepless
Dimensions (W x D x H)	200 x 250 x 90 mm
Weight(kg)	2.9 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	35 W