



LABINCO L27 Large Rotary Mixer

Beschreibung

Large Rotary-mixer Model L27

- Same as model L26 but 50 mm higher and with timer and stronger motor Weight: 8 kg
- Provides rotation action for up to 310 mm diameter accessories
- Large Rotary Mixer for large serial rotation for Bottle / Hybridisation
- With timer and strong motor
- Speed 1-50 rpm, analog
- Provided with stainless steel drip tray
- Quick coupling for exchange of rotating drum or bar
- Dimensions (WxDxH) 540 x 160 x 235 mm
- Weight 8 kg
- Art.nr.27000 (230V) 27010 (115V)
- See also the other models for more options on www.labinco-bv.com

Zustätzliche Informationen

Number of positions	1 bar, drum, bottle holder
Max Quantity	4 kg
Motor Rating input(W)	85 W
Motor Rating output(W)	55 W
Speed display	analog
Speed range(rpm)	1 - 50 rpm

Stirring bar length max(mm)	digital 1s. up to 99h 59m 59s
Speed control	stepless
Dimensions (W x D x H)	540 x 160 x 235 mm
Weight(kg)	8 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	'-220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	100 W