




## LABINCO L28 Test-Tube Rotator

### Beschreibung

#### Model ?? L28

- For mixing, diffusion, dialysing, extraction etc.
  - SPEED 10 rpm With strong, geared motor
  - Angle of rotation can be set between 90 and 27
  - Tubes may be exchanged while rotator is running
  - Comes with line-cord, but without disks
- 
- Test-Tube Rotator for serial rotation of test-tubes 4-40 mm
  - Speed 10 rpm not adjustable
  - Angle of rotation can be set between 90 and 27
  - Dimensions (WxDxH) 200x 200 x 220 mm
  - Weight 2.6 kg
  - Art.nr.28000 (230V) 28010 (115V)
  - See also the other models for more options on [www.labinco-bv.com](http://www.labinco-bv.com)

### Zustätzliche Informationen

|                               |                   |
|-------------------------------|-------------------|
| <b>Number of positions</b>    | 1 for max 2 disks |
| <b>Max Quantity</b>           | 1kg               |
| <b>Motor Rating input(W)</b>  | 30W               |
| <b>Motor Rating output(W)</b> | 25W               |
| <b>Speed range(rpm)</b>       | 10rpm             |
| <b>Speed control</b>          | not adjustable    |

|                                                   |                           |
|---------------------------------------------------|---------------------------|
| <b>Dimensions (W x D x H)</b>                     | 200 x 200 x 220 mm        |
| <b>Weight(kg)</b>                                 | 2.6 kg                    |
| <b>Permissible ambient temperature</b>            | 5 - 60 C                  |
| <b>Permissible relative humidity</b>              | 80 %                      |
| <b>Protection class according to DIN EN 60529</b> | IP 42                     |
| <b>Voltage</b>                                    | 220 - 230 V (115 / 100 V) |
| <b>Frequency</b>                                  | 50/60 Hz                  |
| <b>Power input</b>                                | 35 W                      |