




LABINCO L29 Large Test-Tube Rotator Variable speed

Beschreibung

Model ??? L29

- For mixing, diffusion, dialysing, extraction etc.
 - variable speed 5 to 50 rpm
 - Angle of rotation can be set between 90 and 27
 - Tubes may be exchanged while rotator is running
 - Comes with line-cord, but without disks
-
- Large Test-Tube Rotator for serial rotation of test-tubes 4-40 mm
 - Speed 1-50 rpm analog
 - Angle of rotation can be set between 90 and 27
 - Dimensions (WxDxH) 200x 200 x 270 mm
 - Weight 3.6 kg
 - Art.nr.29000 (100 ??? 230V)

Zustzliche Informationen

Number of positions	1 for max 3 disks
Max Quantity	2 kg
Motor Rating input(W)	40W
Motor Rating output(W)	35 W
Speed display	analog
Speed range(rpm)	1-50 rpm
Speed control	analog

Dimensions (W x D x H)	200 x 200 x 270 mm
Weight(kg)	3.6 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	100 - 230 V)
Frequency	50/60 Hz
Power input	60W