



LABINCO L295 Solubility Index Mixer

Beschreibung

Model L-295

- Solubility Index Mixer is an precision mixer for use primarily in determining the solubility of dry milk powder.
- Specified by the American Dry-Milk Institute in standard method no: B 911 and no: B 913 (and DLG-Prfungsbestimmungen fr Milch and Milch-Erzeugnisse 1965, page 43).
- Strong motor and digital timer
- Special glass jar included (Art.nr.29600)
- Solubility Index Mixer is an precision mixer for use primarily in determining the solubility of dry milk powder.
- Specified by the American Dry-Milk Institute in standard method no: B 911 and no: B 913 (and DLG-Prfungsbestimmungen fr Milch and Milch-Erzeugnisse 1965, page 43).
- Strong motor and digital timer
- Dimensions (WxDxH) 250 x 260 x 380 mm
- Weight 11.2 kg
- Art.nr.29500 (230V) 29510 (115V)
- Special glass jar included (Art.nr.29600)

Zustzliche Informationen

Number of positions	1
Max Quantity	special glass jar
Motor Rating input(W)	60W
Motor Rating output(W)	55 W

Speed range(rpm)	not adjustable
Stirring bar length max(mm)	digital 1s. up to 99h 59m 59s
Speed control	not adjustable
Adjustable safety circuit min	2 x
Dimensions (W x D x H)	250 x 260 x 380 mm
Weight(kg)	11.2 kg
Permissible ambient temperature	5 - 60 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	65 W