



## LABINCO L295 Solubility Index Mixer

### Beschreibung

#### Model L-295

- Solubility Index Mixer is an precision mixer for use primarily in determining the solubility of dry milk powder.
- Specified by the American Dry-Milk Institute in standard method no: B 911 and no: B 913 (and DLG-Prfungsbestimmungen fr Milch and Milch-Erzeugnisse 1965, page 43).
- Strong motor and digital timer
- Special glass jar included ( Art.nr.29600 )
- Solubility Index Mixer is an precision mixer for use primarily in determining the solubility of dry milk powder.
- Specified by the American Dry-Milk Institute in standard method no: B 911 and no: B 913 (and DLG-Prfungsbestimmungen fr Milch and Milch-Erzeugnisse 1965, page 43).
- Strong motor and digital timer
- Dimensions (WxDxH) 250 x 260 x 380 mm
- Weight 11.2 kg
- Art.nr.29500 (230V) 29510 (115V)
- Special glass jar included ( Art.nr.29600 )

### Zustzliche Informationen

|                        |                   |
|------------------------|-------------------|
| Number of positions    | 1                 |
| Max Quantity           | special glass jar |
| Motor Rating input(W)  | 60W               |
| Motor Rating output(W) | 55 W              |

|                                                   |                               |
|---------------------------------------------------|-------------------------------|
| <b>Speed range(rpm)</b>                           | not adjustable                |
| <b>Stirring bar length max(mm)</b>                | digital 1s. up to 99h 59m 59s |
| <b>Speed control</b>                              | not adjustable                |
| <b>Adjustable safety circuit min</b>              | 2 x                           |
| <b>Dimensions (W x D x H)</b>                     | 250 x 260 x 380 mm            |
| <b>Weight(kg)</b>                                 | 11.2 kg                       |
| <b>Permissible ambient temperature</b>            | 5 - 60 C                      |
| <b>Permissible relative humidity</b>              | 80 %                          |
| <b>Protection class according to DIN EN 60529</b> | IP 42                         |
| <b>Voltage</b>                                    | 220 - 230 V (115 / 100 V)     |
| <b>Frequency</b>                                  | 50/60 Hz                      |
| <b>Power input</b>                                | 65 W                          |