



LABINCO L30 Hotplate 20 ltr. 325°C

Beschreibung

Model L-30

- Top plate: Aluminum 195 mm diam.
- With Stainless Steel efficient reflectors
- Thermostatic controlled 50-320C
- Housing: 2 piece casted aluminum
- Electrostatically applied powder coated

Hotplate 195 mm diam. Aluminium

- Hotplate temperature 50-325C
- Dimensions (WxDxH) 200 x 250 x 135 mm
- Weight 3,8 kg
- Art.nr.30000 (230V) 30010 (115V)

Zustzliche Informationen

Number of positions	1
Max Quantity	20 ltr
Heat output	600 W
Heating rate	7 K/min
Heating temperature range	room temp. - 325 C
Heat control	analog
Heat control accuracy	1 K
Temperature constancy in medium	1 K

Dimensions (W x D x H)	200 x 250 x 135 mm
Weight(kg)	3.8 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	600 W