



LABINCO L36 Flask Heater 450°C

Beschreibung

Model L-36

- Heatingplate: 140 mm (steel)
- Controlled 50 °C - 450°C
- Heating: 450 W
- Hotplate 85 mm diam. steel coated
- Hotplate temperature 50-450°C
- Dimensions (WxDxH) 172 x 286 x 125 mm
- Weight 1.8 kg
- Art.nr.36000 (230V)
- Aluminium flask holders

Zusätzliche Informationen

Number of positions	1
Max Quantity	100-1000 ml
Heat output	450 W
Heating rate	7 K/min
Heating temperature range	room temp. - 450 C
Heat control	analog
Heat control accuracy	1 K
Temperature constancy in medium	1 K
Dimensions (W x D x H)	172 x 286 x 125 mm

Weight(kg)	1.8kg
Permissible ambient temperature	5 - 60 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	450 W