



LABINCO L364 4-Pos. Flask Heater 450Â°C

Beschreibung

- 4 x hotplate 85 mm diam. steel coated
- Hotplate temperature 50-450C
- Dimensions (WxDxH) 600 x 240 x 130 mm
- Weight 8.2 kg
- Art.nr.36400 (230V)
- Aluminium flask holders

- 4 x hotplate 85 mm diam. steel coated
- Hotplate temperature 50-450C
- Dimensions (WxDxH) 600 x 240 x 130 mm
- Weight 8.2 kg
- Art.nr.36400 (230V)
- Aluminium flask holders

Zustzliche Informationen

Number of positions	4
Max Quantity	100-1000 ml (see flask holders)
Heat output	450 W (1 Position)
Heating rate	7 K/min
Heating temperature range	room temp. - 450 C
Heat control	analog
Heat control accuracy	1 K
Temperatureconstancyinmedium	1 K

Dimensions (W x D x H)	600 x 240 x 130 mm
Weight(kg)	8.2 kg
Permissible ambient temperature	5 - 60 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	1800 W