



## LABINCO L366 6-Pos. Flask Heater 450°C

### Beschreibung

- 6 x hotplate 85 mm diam. steel coated
- Hotplate temperature 50-450C
- Dimensions (WxDxH) 900 x 240 x 130 mm
- Weight 11.2 kg
- Art.nr.36600 (230V)
- Aluminium flask holders

- 6 x hotplate 85 mm diam. steel coated
- Hotplate temperature 50-450C
- Dimensions (WxDxH) 900 x 240 x 130 mm
- Weight 11.2 kg
- Art.nr.36600 (230V)
- Aluminium flask holders

### Zustzliche Informationen

|                              |                                 |
|------------------------------|---------------------------------|
| Number of positions          | 6                               |
| Max Quantity                 | 100-1000 ml (see flask holders) |
| Heat output                  | 450 W (1 Position )             |
| Heating rate                 | 7 K/min                         |
| Heating temperature range    | room temp. - 450 C              |
| Heat control                 | analog                          |
| Heat control accuracy        | 1 K                             |
| Temperatureconstancyinmedium | 1 K                             |

|                                            |                    |
|--------------------------------------------|--------------------|
| Dimensions (W x D x H)                     | 900 x 240 x 130 mm |
| Weight(kg)                                 | 11.2 kg            |
| Permissible ambient temperature            | 5 - 60 C           |
| Permissible relative humidity              | 80 %               |
| Protection class according to DIN EN 60529 | IP 42              |
| Voltage                                    | 220 - 230 V        |
| Frequency                                  | 50/60 Hz           |
| Power input                                | 2700 W             |