



LABINCO L366 6-Pos. Flask Heater 450°C

Beschreibung

- 6 x hotplate 85 mm diam. steel coated
 - Hotplate temperature 50-450C
 - Dimensions (WxDxH) 900 x 240 x 130 mm
 - Weight 11.2 kg
 - Art.nr.36600 (230V)
 - Aluminium flask holders
-
- 6 x hotplate 85 mm diam. steel coated
 - Hotplate temperature 50-450C
 - Dimensions (WxDxH) 900 x 240 x 130 mm
 - Weight 11.2 kg
 - Art.nr.36600 (230V)
 - Aluminium flask holders

Zustzliche Informationen

Number of positions	6
Max Quantity	100-1000 ml (see flask holders)
Heat output	450 W (1 Position)
Heating rate	7 K/min
Heating temperature range	room temp. - 450 C
Heat control	analog
Heat control accuracy	1 K
Temperatureconstancyinmedium	1 K

Dimensions (W x D x H)	900 x 240 x 130 mm
Weight(kg)	11.2 kg
Permissible ambient temperature	5 - 60 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	2700 W