




LABINCO L51 12-channel fixed speed pump

Beschreibung

- Pump for up to 12 pump tubes 0.03-3.70 ml/min
 - 1 pump bed with continuously adjust. tubing tension
 - With 6 rollers, minizing pulsation and maximizing back-pressure level
 - Speed 10 rpm not adjustable
 - Dimensions (WxDxH) 260x 130 x 120 mm
 - Weight 7 kg
 - Art.nr.51000 (230V) 51800 (115V)
-
- Pump for up to 12 pump tubes 0.03-3.70 ml/min
 - 1 pump bed with continuously adjust. tubing tension
 - With 6 rollers, minizing pulsation and maximizing back-pressure level
 - Speed 10 rpm not adjustable
 - Dimensions (WxDxH) 260x 130 x 120 mm
 - Weight 7 kg
 - Art.nr.51000 (230V) 51800 (115V)

Zustzliche Informationen

Number of positions	12
Max Quantity	0.03-3.70 ml/min
Motor Rating input(W)	75W
Motor Rating output(W)	55 W
Speed range(rpm)	10 rpm not adjustable
Speed control	not adjustable

Dimensions (W x D x H)	260 x 130 x 120 mm
Weight(kg)	7 kg
Permissible ambient temperature	5 - 60 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	75 W