



LABINCO L71 Basic Magnetic Stirrer 10 ltr.

Beschreibung

- Topplate 155 mm diam. Stainless steel
- Stirring capacity 10 Ltr
- Dimensions (WxDxH) 172 x 286 x 80 mm
- Weight 2.2 kg
- Art.nr.71000 (230V) 71010 (115V)
- Incl. spinbar

- Topplate 155 mm diam. Stainless steel
- Stirring capacity 10 Ltr
- Dimensions (WxDxH) 172 x 286 x 80 mm
- Weight 2.2 kg
- Art.nr.71000 (230V) 71010 (115V)
- Incl. spinbar

Zustzliche Informationen

Number of positions	1
Max Quantity	10 ltr
Motor Rating input(W)	35 W
Motor Rating output(W)	15 W
Speed display	analog
Speed range(rpm)	50-2500rpm
Stirring bar length max(mm)	80 mm
Speed control	stepless

Dimensions (W x D x H)	172 x 286 x 80 mm
Weight(kg)	2.2 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	35 W