



LABINCO L73 Large Magnetic Stirrer 20 ltr.

Beschreibung

- Topplate 155 mm diam. Stainless steel
- Stirring capacity 20 Ltr
- Dimensions (WxDxH) 172 x 286 x 80 mm
- Weight 2.4 kg
- Art.nr.73000 (230V)
- Incl. spinbar
- Topplate 155 mm diam. Stainless steel
- Stirring capacity 20 Ltr
- Dimensions (WxDxH) 172 x 286 x 80 mm
- Weight 2.4 kg
- Art.nr.73000 (230V)
- Incl. spinbar

Zustzliche Informationen

Number of positions	1
Max Quantity	20 ltr
Motor Rating input(W)	50W
Motor Rating output(W)	20 W
Speed display	analog
Speed range(rpm)	50 - 1250 rpm
Stirring bar length max(mm)	80 mm
Speed control	stepless

Dimensions (W x D x H)	172 x 286 x 80 mm
Weight(kg)	2.4 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	35 W