



LABINCO L746 6 position Magnetic Stirrer 6 x 5 ltr.

Beschreibung

- Topplate 495 x 360 mm diam. Stainles steel 50-1250 rpm
 - 6 positions in two rows of three.
 - Stirring capacity 6 x 5 Ltr
 - Dimensions (WxDxH) 510 x 405 x 145 mm
 - Weight 12 kg
 - Art.nr.74600 (230V)
 - Incl. 6 x spinbar
-
- Topplate 495 x 360 mm diam. Stainles steel 50-1250 rpm
 - 6 positions in two rows of three.
 - Stirring capacity 6 x 5 Ltr
 - Dimensions (WxDxH) 510 x 405 x 145 mm
 - Weight 12 kg
 - Art.nr.74600 (230V)
 - Incl. 6 x spinbar

Zustzliche Informationen

Number of positions	6 positions in two rows of three.
Max Quantity	6 x 5 ltr
Motor Rating input(W)	6 x 35W
Motor Rating output(W)	6 x 15 W
Speed display	analog
Speed range(rpm)	50-2500rpm

Stirring bar length max(mm)	60 mm
Speed control	stepless
Dimensions (W x D x H)	510 x 405 x 145 mm
Weight(kg)	12 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	220 W