



LABINCO L80 Hotplate 10 ltr. 325°C

Beschreibung

Model L-80 St.steel

- Hotplate: Aluminium 145 mm diam.
- With Stainless Steel efficient reflectors
- Thermostatic controlled 50-320C
- Housing from casted aluminum
- Electrostatically applied powder coated
- M-10 threaded support rod connection
- Hotplate 145 mm diam. Aluminium
- Hotplate temperature 50-325C
- Dimensions (WxDxH) 172 x 286 x 125 mm
- Weight 2 kg
- Art.nr.80000 (230V) 80010 (115V)

Zustzliche Informationen

Number of positions	1
Max Quantity	10 ltr
Heat output	600 W
Heating rate	7 K/min
Heating temperature range	room temp. - 325 C
Heat control	analog
Heat control accuracy	1 K

Temperatureconstancyinmedium	1 K
Dimensions (W x D x H)	172 x 286 x 125 mm
Weight(kg)	2 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	600 W