



## LABINCO LD-48 Shaker

### Beschreibung

Labinco LD-48 Laboratory Shaker.Labinco LD-48 Laboratory Shaker.

### Zustzliche Informationen

|                                                   |                                                  |
|---------------------------------------------------|--------------------------------------------------|
| <b>Number of positions</b>                        | See LD-48 parts (for a lot of different plates). |
| <b>Motor Rating input(W)</b>                      | 90 W                                             |
| <b>Motor Rating output(W)</b>                     | 60 W                                             |
| <b>Speed display</b>                              | Digital                                          |
| <b>Speed range(rpm)</b>                           | 75 - 1250 rpm                                    |
| <b>Speed control</b>                              | Digital                                          |
| <b>Dimensions (W x D x H)</b>                     | 310 x 340 x 140 mm                               |
| <b>Weight(kg)</b>                                 | 13 kg                                            |
| <b>Permissible ambient temperature</b>            | 5 - 40 C                                         |
| <b>Permissible relative humidity</b>              | 80%                                              |
| <b>Protection class according to DIN EN 60529</b> | IP 42                                            |
| <b>Voltage</b>                                    | 100 - 230 V                                      |
| <b>Frequency</b>                                  | 50/60 Hz                                         |
| <b>Power input</b>                                | 125 W                                            |