




LABINCO LD14 Large Magnetic Stirrer 60 ltr.

Beschreibung

Large Stirrer 35 liter LD-14

- Optical speed-detection
- Reversible stirring
- Topplate st.st.
- Supplied with a teflon magnetic stirring bar line cord and manual
- Topplate 320 x 340 mm Stainless steel 30-1000 rpm digital
- Stirring capacity 60 Ltr
- Dimensions (WxDxH) 320 x 340 x 105 mm
- Weight 5.2 kg
- Art.nr.14000 (230V)
- Incl. spinbar

Zustzliche Informationen

Number of positions	1
Max Quantity	60 ltr
Motor Rating input(W)	50 W
Motor Rating output(W)	20 W
Speed display	digital
Speed range(rpm)	30 - 1000rpm
Stirring bar length max(mm)	120 mm
Speed control	stepless

Dimensions (W x D x H)	320 x 340 x 105 mm
Weight(kg)	5.2 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	50 W