




LABINCO LD30 Hotplate 20 ltr. 350°C

Beschreibung

- Hotplate 195 mm diam. Aluminium
- Hotplate temperature 50-350C
- Dimensions (WxDxH) 200 x 250 x 135 mm
- Weight 3.8 kg
- Art.nr.30001 (230V)
- Incl. temperature sensor

- Hotplate 195 mm diam. Aluminium
- Hotplate temperature 50-350C
- Dimensions (WxDxH) 200 x 250 x 135 mm
- Weight 3.8 kg
- Art.nr.30001 (230V)
- Incl. temperature sensor

Zustzliche Informationen

Number of positions	1
Max Quantity	20 ltr
Heat output	600 W
Heating rate	7 K/min
Heating temperature range	room temp. - 350 C
Heat control	digital
Heat control accuracy	1 K
Connection for ext. temperature sensor	Termocouple

Control accuracy with sensor	1 K
Temperatureconstancyinmedium	1 K
Dimensions (W x D x H)	200 x 250 x 135 mm
Weight(kg)	3.8 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	610 W