




LABINCO LD31 Large Hotplate 350°C

Beschreibung

- Hotplate 460 x 350 mm Aluminium
- Hotplate temperature 30-350C digital
- Dimensions (WxDxH) 470 x 370 x 145 mm
- Weight 9.1 kg
- Temperature sensor on request
- Art.nr.31001 (230V)

- Hotplate 460 x 350 mm Aluminium
- Hotplate temperature 30-350C digital
- Dimensions (WxDxH) 470 x 370 x 145 mm
- Weight 9.1 kg
- Temperature sensor on request
- Art.nr.31001 (230V)

Zustzliche Informationen

Number of positions	1 large hotplate for more pos.
Max Quantity	1 ml -50 ltr or more
Heat output	2000 W
Heating rate	7 K/min
Heating temperature range	room temp. - 350 C
Heat control	digital
Heat control accuracy	1 K
Connection for ext. temperature sensor	on request

Control accuracy with sensor	1 K
Temperatureconstancyinmedium	1 K
Dimensions (W x D x H)	470 x 370 x 145 mm
Weight(kg)	9.1 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	2000 W