




LABINCO LD39 Digital Hotplate Magnetic Stirrer 60 ltr. 450Â°C

Beschreibung

- Hotplate 300 mm diam. Aluminium 10 â?? 1000 rpm
 - Hotplate temperature 25-450C (digital)
 - Stirring capacity 60 Ltr
 - Dimensions (WxDxH) 320 x 340 x 155 mm
 - Weight 9.1 kg
 - Art.nr.39000 (230V)
 - Temperature sensor on request
 - Incl. spinbar
-
- Hotplate 300 mm diam. Aluminium 10 â?? 1000 rpm
 - Hotplate temperature 25-450C (digital)
 - Stirring capacity 60 Ltr
 - Dimensions (WxDxH) 320 x 340 x 155 mm
 - Weight 9.1 kg
 - Art.nr.39000 (230V)
 - Temperature sensor on request
 - Incl. spinbar

Zustzliche Informationen

Number of positions	1
Max Quantity	60 ltr
Motor Rating input(W)	100 W
Motor Rating output(W)	45 W
Speed display	analog

Speed range(rpm)	10 - 1000 rpm
Stirring bar length max(mm)	120 mm
Heat output	2000 W
Heating rate	7 K/min
Heating temperature range	room temp. - 450 C
Heat control	digital
Heat control accuracy	1 K
Speed control	stepless
Connection for ext. temperature sensor	on request
Control accuracy with sensor	1 K
Temperature constancy in medium	1 K
Dimensions (W x D x H)	320 x 340 x 155 mm
Weight(kg)	9.1 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
RS232 interface	on request
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	2100 W