




LABINCO LD76 Digital Rotary Mixer

Beschreibung

- Digital Rotary Mixer for small serial rotation of test-tubes 4-27 mm
 - With timer
 - Memory function and Re-start function
 - Speed 1-99 rpm, digital
 - Quick coupling for exchange of rotating drum or bar
 - Dimensions (WxDxH) 490 x 160 x 205 mm
 - Weight 4.5 kg
 - Art.nr.76000 (230V) 76010 (115V)
-
- Digital Rotary Mixer for small serial rotation of test-tubes 4-27 mm
 - With timer
 - Memory function and Re-start function
 - Speed 1-99 rpm, digital
 - Quick coupling for exchange of rotating drum or bar
 - Dimensions (WxDxH) 490 x 160 x 205 mm
 - Weight 4.5 kg
 - Art.nr.76000 (230V) 76010 (115V)

Zusätzliche Informationen

Number of positions	1 bar or drum
Max Quantity	1 kg
Motor Rating input(W)	50 W
Motor Rating output(W)	35 W
Speed display	digital

Speed range(rpm)	1 - 99 rpm
Stirring bar length max(mm)	digital 1s. up to 99h 59m 59s
Speed control	digital
Dimensions (W x D x H)	540 x 160 x 235 mm
Weight(kg)	8 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	60 W