




LABINCO LD79 Digital Test-Tube Rotator

Beschreibung

- Digital Test-Tube Rotator for serial rotation of test-tubes 4-40 mm
 - For rotation disks see www.labinco-bv.com
 - With timer
 - Memory function and Re-start function
 - Speed 1-99 rpm, digital
 - Angle of rotation can be set between 90 and 27
 - Dimensions (WxDxH) 200x 200 x 250 mm
 - Weight 3.5 kg
 - Art.nr.79000 (230V) 79010 (115V)
-
- Digital Test-Tube Rotator for serial rotation of test-tubes 4-40 mm
 - For rotation disks see www.labinco-bv.com
 - With timer
 - Memory function and Re-start function
 - Speed 1-99 rpm, digital
 - Angle of rotation can be set between 90 and 27
 - Dimensions (WxDxH) 200x 200 x 250 mm
 - Weight 3.5 kg
 - Art.nr.79000 (230V) 79010 (115V)

Zustzliche Informationen

Number of positions	1 for max 2 disks
Max Quantity	0.5 kg
Motor Rating input(W)	50 W

Motor Rating output(W)	35 W
Speed display	digital
Speed range(rpm)	1 - 99 rpm
Stirring bar length max(mm)	digital 1s. up to 99h 59m 59s
Speed control	digital
Dimensions (W x D x H)	200 x 200 x 250 mm
Weight(kg)	3.5 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V (115 / 100 V)
Frequency	50/60 Hz
Power input	60 W