



LABINCO LD80-S Digital contact Thermometer Stainless steel probe

Beschreibung

- Waterproof precision regulating thermometer (pt100) for temperatures between -50C to +300C
 - Overshooting of the temperature is minimized by a fuzzy control
 - Digitally adjustable
 - Stainless steel probe for use with especially aggressive substances
 - Support rod with clamps for probe (not included) Art.nr.80900
 - 3 mm and 210 mm long (300 mm with handle)
 - Dimensions (WxDxH) 55 x 25 x 145 mm
 - Weight 0.4 kg
 - Art.nr.80000 (9â?¡12 Vdc)
-
- Waterproof precision regulating thermometer (pt100) for temperatures between -50C to +300C
 - Overshooting of the temperature is minimized by a fuzzy control
 - Digitally adjustable
 - Stainless steel probe for use with especially aggressive substances
 - Support rod with clamps for probe (not included) Art.nr.80900
 - 3 mm and 210 mm long (300 mm with handle)
 - Dimensions (WxDxH) 55 x 25 x 145 mm
 - Weight 0.4 kg
 - Art.nr.80000 (9â?¡12 Vdc)

Zustzliche Informationen

Number of positions	1
Speed display	digital adjustable
Heat control	-50C to +300C

Dimensions (W x D x H)	55 x 25 x 145 mm
Weight(kg)	0.4 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80 %
Protection class according to DIN EN 60529	IP 42
Voltage	9...12 Vdc
Power input	10 W