



LABINCO LD80-S Digital contact Thermometer Stainless steel probe

Beschreibung

- Waterproof precision regulating thermometer (pt100) for temperatures between -50C to +300C
- Overshooting of the temperature is minimized by a fuzzy control
- Digitally adjustable
- Stainless steel probe for use with especially aggressive substances
- Support rod with clamps for probe (not included) Art.nr.80900
- 3 mm and 210 mm long (300 mm with handle)
- Dimensions (WxDxH) 55 x 25 x 145 mm
- Weight 0.4 kg
- Art.nr.80000 (9â?l12 Vdc)
- Waterproof precision regulating thermometer (pt100) for temperatures between -50C to +300C
- Overshooting of the temperature is minimized by a fuzzy control
- Digitally adjustable
- Stainless steel probe for use with especially aggressive substances
- Support rod with clamps for probe (not included) Art.nr.80900
- 3 mm and 210 mm long (300 mm with handle)
- Dimensions (WxDxH) 55 x 25 x 145 mm
- Weight 0.4 kg
- Art.nr.80000 (9â?l12 Vdc)

Zustzliche Informationen

Anzahl der Positionen	1
Geschwindigkeitsanzeige	digital adjustable
Wrmesteuerung	-50C to +300C

Mae (L x B x H)	55 x 25 x 145 mm
Gewicht	0.4 kg
Zulssige Umgebungstemperatur	5 - 40 C
Max. Luftfeuchtigkeit	80 %
Sicherheitsklasse DIN EN 60529	IP 42
Spannung	9...12 Vdc
Leistungsaufnahme	10 W