



LABINCO LD816 with probe

Beschreibung

Zustzliche Informationen

Number of positions	4 or 6 Bench of Square models
Max Quantity	4 or 6 x 5 ltr
Motor Rating input(W)	4 or 6 x 35 W
Motor Rating output(W)	4 or 6 x 15 W
Speed display	Analog
Speed range(rpm)	50 - 1250 rpm
Stirring bar length max(mm)	60mm
Heat output	4 or 6 x 300 W
Heating rate	7 K/min
Heating temperature range	Room temp up to 350 C
Heat control	Digital
Heat control accuracy	1 K
Speed control	Stepless
Adjustable safety circuit min	0 C
Adjustable safety circuit max	350 C
Connection for ext. temperature sensor	On request
Control accuracy with sensor	1 K
Temperatureconstancyinmedium	1 K
Dimensions (W x D x H)	835 x 245 140 mm (6 pos.)

Weight(kg)	14 kg
Permissible ambient temperature	4 -50 C
Permissible relative humidity	80%
Protection class according to DIN EN 60529	IP42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	6 pox. 2050 W (4 pos. 1350 W)