



LABINCO Magnetic stirrer with stainless steel heating

Beschreibung

This heating with stirrer system is used to mix the samples in a glass jars in a automated robotic process to quickly heat and mix glass jars with samples. It is custom made and developed for this customer, completely made out of stainless steel with drip tray.

All our custom made products are designed / developed and created in our own factory in The Netherlands.

Custom made product, let us know your wishes.

Zusätzliche Informationen

| | |
|--------------------------------|------------|
| Anzahl der Positionen | 4 |
| Drehzahlbereich (U/min) | 100 - 3000 |
| Spannung | 230 V |
| Leistungsaufnahme | 2500 W |