




LABINCO Special big rolling bench

Beschreibung

Zustzliche Informationen

Number of positions	1 up to 8
Motor Rating input(W)	100 W
Motor Rating output(W)	75 W
Speed display	Analog
Speed range(rpm)	5 - 240 rpm
Speed control	Stepless
Dimensions (W x D x H)	950 x 550 x 135
Weight(kg)	15.5
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80%
Protection class according to DIN EN 60529	IP 42
Voltage	220 - 230 V
Frequency	50/60 Hz
Power input	150 W