




LABINCO very strong Power Mixer

Beschreibung

Zustzliche Informationen

Number of positions	1
Motor Rating input(W)	90 W
Motor Rating output(W)	60 W
Speed display	Analog
Speed range(rpm)	50 - 2500 rpm
Stirring bar length max(mm)	45 mm or more (on request)
Speed control	Analog (stepless)
Dimensions (W x D x H)	200 x 270 x 170 mm
Weight(kg)	14 kg
Permissible ambient temperature	5 - 40 C
Permissible relative humidity	80%
Protection class according to DIN EN 60529	IP 42
Voltage	100 - 230 V
Frequency	50/60 Hz
Power input	125 W